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Arrow® T-POD™ Pelvic Stabilization Device

Rx only:

Indications for Use:

Provides circumferential compression to the pelvis in patients with suspected pelvic fracture for pelvic stabilization which may reduce blood loss and pain.

Disposal:


Where applicable dispose of according to local regulations for biohazardous waste.

When Used In X-Ray:

Detachable X-Ray Detectable Label is included to identify that T-POD is in place on patient.

⚠ General Warnings

1. Warning: Single use: Do not reuse, reprocess or sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Photos are for illustration purposes.

Instruction For Use:

1. Slide the Belt under the patient and into position under the pelvis. It is recommended to remove patient clothing prior to applying T-POD.
2. Trim the Belt leaving a 6-8" (15-20 cm) gap over the center of the abdomen. Note: If unable to trim, the Belt can be folded under.
3. Apply the Velcro®-backed Pulley System to the Belt on each side of the gap.
4. Draw the Pull Tab, creating simultaneous circumferential compression.
5. Wrap the cord around the Pulley System hooks and secure the Velcro®-backed Pull Tab to the Pulley System fabric.
6. Record the date and time of application.

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Teleflex, das Teleflex-Logo, Arrow, das Arrow-Logo und T-POD sind Marken oder eingetragene Marken von Teleflex Incorporated oder verbundenen Unternehmen in den USA und/oder anderen Ländern. Andere Namen sind eventuell Marken des jeweiligen Eigentümers. © 2019 Teleflex Incorporated. Alle Rechte vorbehalten.