

Larynx-Masken

- Entlüften
- Gel
- Einführen bis Halter 1,5-2cm zur Lippe
- Fixieren
- Blocken (ca. 60cm H₂O)
- Bubble Test, Abhören, Kapnometrie
- Magensonde
- Larynxtubus !!!
- https://www.youtube.com/watch?v=7WAu_4-47Lw
- lma-supreme.pdf

A

From:

<https://www.rdbuch.lima-city.de/> - Wiki

Permanent link:

<https://www.rdbuch.lima-city.de/doku.php?id=praxisanleiter:einweisungen:larynx-masken>

Last update: **2024/06/03 15:58**

